



ECHELON[◇] Primary Hip System

Surgical Technique

Harry B. Skinner, MD, PhD
Professor and Chairman
Department of Orthopaedic Surgery
University of California, Irvine
Orange, California

James P. Waddell, MD, FRCS(C)
A.J. Latner, Professor & Chairman
Division of Orthopaedic Surgery
University of Toronto
Toronto, Ontario, Canada


Douglas A. Becker, MD
Assistant Professor
University of Minnesota
School of Medicine –
Department of Orthopaedics
Minneapolis, Minnesota

Paul S. Lux, MD
The Orthopedic Center of St. Louis
Barnes-Jewish West County Hospital
St. Louis, Missouri


Nota Bene

The technique description herein is made available to the healthcare professional to illustrate the authors' suggested treatment for the uncomplicated procedure. In the final analysis, the preferred treatment is that which addresses the needs of the patient.


Basic steps



Ream to size



Broach and trial



Implant stem

Porous femoral implants

Neck geometry –

Circulotrapezoidal neck provides increased range of motion compared to a circular neck of the same strength. Polished surface improves fatigue strength.

Driving platform –

The ECHELON® implants feature a threaded driving platform with an elliptical slot for rotational and axial implant control during insertion.

Lateral proximal flare –

ECHELON has a 3° proximal anterior/posterior flare to improve proximal fill, without preventing implant seating.


Optimized 12/14 taper –

Modular femoral heads are locked onto the stem with an optimized 12/14 taper which contributes to an increased range of motion.

Distal bullet tip –

The bullet tip reduces the stress between the distal implant tip and the bone to minimize end-of-stem thigh pain.

Standard offset



Shoulder relief –

The lateral shoulder is rounded to minimize the risk of fracturing the greater trochanter during stem insertion.

Neck offset options –

Standard and high offset options are available to ensure the appropriate joint tension.

Porous coating –

ROUGHCOAT® porous coating increases the friction between the implant and bone, improving implant stability and providing a porous surface for bone ingrowth.

Material –

All ECHELON® implants are manufactured from Cobalt Chromium allowing for extensive porous coating of the stem.

Size range –


ECHELON Primary stems are offered in 1mm increments to minimize bone removal and provide optimum canal fill.

Distal slot –

The distal slot eases stem insertion, reduces the risk of fracture and reduces distal stem stiffness.

Distal flutes –

The ECHELON system offers distal flutes to increase rotational stability.



High offset

Implant specifications

General specifications

Cobalt chromium material

Neck shaft angle 131°

Standard collar shaft angle 50°

Primary stem length 130-160mm*

Porous-coating length 90-108mm**

The broach is 0.5mm smaller than the implant.

Distal flute diameter is 0.25mm larger than porous coated cylindrical diameter.

* Stem length is measured from the collar to the distal tip

** Porous coating length is measured from the shoulder to the distal end of the coating

	Standard offset	Neck length (mm)												
		Size	-3		+0		+4		+8		+12		+16	
High offset		11	24	—	27	—	31	—	35	—	39	—	43	—
		12	24	28	27	31	31	35	35	39	39	43	43	47
		13-14	27	33	30	36	34	40	38	44	42	48	46	52
		15-17	31	36	34	39	38	43	42	47	46	51	50	55
		18-19*	34	39	37	42	41	46	45	50	49	54	53	58

Neck offset (mm)													
Size	-3		+0		+4		+8		+12		+16		
11	32	—	34	—	37	—	40	—	43	—	46	—	
12	32	38	34	40	37	43	40	46	43	49	46	52	
13-14	35	43	37	45	40	48	43	51	46	54	49	57	
15-17	38	46	40	48	43	51	46	54	49	57	52	60	
18-19*	41	49	43	51	46	54	49	57	52	60	55	63	


Neck height (mm)													
Size	-3		+0		+4		+8		+12		+16		
11	25	—	27	—	30	—	32	—	35	—	37	—	
12	25	25	27	27	30	30	32	32	35	35	37	37	
13-14	28	28	30	30	33	33	35	35	38	38	40	40	
15-17	30	30	32	32	35	35	37	37	40	40	42	42	
18-19*	34	34	36	36	39	39	41	41	44	44	46	46	

* 18-19 available as special request

Length measurements

Standard/high offset size	Stem length	Porous coating length
11-12	130mm	90mm
13-14	140mm	96mm
15-17	150mm	102mm
18-19*	160mm	108mm

*18-19 available as special request




Preoperative planning

Preoperative planning is essential to ensure appropriate implant selection and reaming. It may also be valuable in determining management of leg length discrepancy. Both standard offset and high offset implants are available and templates for these implants should be utilized to determine the optimum offset for each individual patient.

Both an anteroposterior radiograph of the pelvis with the hips in neutral rotation and a lateral hip radiograph optimize preoperative templating. The proximal one-third of the femur should be visible on these radiographs.

In order to determine appropriate offset and management of leg length discrepancy, it is necessary that the femoral templates be utilized in conjunction with acetabular templates appropriate for the implant that has been selected for the acetabular reconstruction.

The femoral templates provided allow for different neck offset options in addition to varying head depths and head diameters to ensure comprehensive selection of implants to deal with variations in femoral and acetabular anatomy.



Femoral neck osteotomy

Step 1

The point of the femoral neck resection should be marked with electrocautery corresponding to both the preoperative templating and the intraoperative measurement. This will confirm that the level of the femoral neck resection is appropriate and will re-establish the desired leg length of the proximal femur.

An osteotomy guide is available for proximal bone resection. Resect the proximal bone by cutting through the angled slot. The osteotomy guide has a vertical scale in 5mm increments to help gauge neck height. Osteotomize the femoral neck (Figures 1A and 1B).




Figure 1A




Figure 1B

Acetabular and femoral canal preparation

Step 2

The acetabulum should be prepared in the recommended fashion for the acetabular component to be utilized according to preoperative planning and templating.




Figure 3A

Step 3

Remove remnants of the femoral neck and open the medullary canal using the box osteotome (Figure 3A). Use the canal finder and modular T-handle for initial femoral reaming (Figure 3B).

Note It is important to stay lateral with both the box osteotome and canal finder. Care should be taken to ensure that the initial reaming track into the femur is in neutral alignment with the femoral axis.

Do not pressurize intramedullary contents with canal finder by inserting too rapidly.

Caution Take care when handling reamers and broaches as they are sharp and may damage surgical gloves and soft tissue.




Figure 3B

Femoral reaming

Step 4

Rigid femoral reamers in 0.5mm increments are available.

The stem size is measured at the maximum diameter of the distal porous coating. The maximum diameter of the flutes is 0.25mm larger than the diameter of the porous coating.

Start reaming with a reamer 4 to 6mm smaller than the templated size or a reamer that has little or no resistance in the femoral canal.

For a line-to-line fit, ream the canal in 0.5mm increments until the last reamer matches the selected implant size. The canal can also be reamed 0.5mm smaller than the size for a tighter distal fit. The final reamer size should be based on bone quality, anatomy and surgeon preference.

Note The flutes on the distal stem are 0.25mm larger than the porous coated diameter. Therefore, reaming line-to-line will produce a 0.25mm press fit in the distal fluted region of the stem.

The stem length is measured from the collar to the distal tip of the implant. Reaming depth is also measured from the collar to the distal tip of the implant (Figure 4A and Figure 4B). Use the implant reaming chart to determine the reaming depth for the porous implants. Seat the reamer to the appropriate depth mark on each reamer.




Figure 4A




Figure 4B

Implant reaming chart

Stem	Ream depth from medial resection level
11-12	130mm
13-14	140mm
15-17	150mm
18-19*	160mm

* 18-19 available as special request

Note ECHELON® reamers have multiple depth marks. The distal-most mark indicates the Primary stem reaming depth. Other depth markings indicate reaming depths for ECHELON revision implants.

Femoral broaching and calcar preparation

Step 5

Attach the broach handle to the femoral broach. Begin broaching two sizes smaller than the size of the last femoral reamer. The broach should be seated to the depth of the medial resection line as shown in Figure 5. The final broach should match the size of the selected implant. The femoral broaches are 0.5mm smaller than the porous coating level of the implant.

Be sure to stay lateral with the smaller broaches to avoid varus broaching.



Figure 5

Step 6

A calcar reamer is available. With the final broach fully seated, remove the broach handle. Place the calcar reamer over the broach post and ream the calcar flush with the top of the broach (Figure 6). This will ensure uniform contact between the collar of the prosthesis and the calcar. **This can be verified by placing the trial neck on the broach.**




Figure 6

Trialing

Step 7

After calcar reaming, place the matching trial neck onto the broach post (as determined by pre-op templating). Fully engage the desired trial femoral head on the trial neck and reduce the hip to assess stability and range of motion (Figure 7).




Figure 7

Femoral neck length options						
Global trial head color	22mm	28mm	32mm	36mm	40mm	44mm
Green	—	XS/-3	XS/-3	XS/-3	XS/-3	XS/-3
Rust	S/+0	S/+0	S/+0	S/+0	S/+0	S/+0
Brown	M/+4	M/+4	M/+4	M/+4	M/+4	M/+4
Gray	L/+8	L/+8	L/+8	L/+8	L/+8	L/+8
Blue	XL/+12*	XL/+12*	XL/+12*	XL/+12	XL/+12	XL/+12
Black	—	XXL/+16*	XXL/+16*	—	—	—

*Skirted femoral head

Trial reduction

Step 8

Reduce the hip and evaluate in the following ways:

- 1 Soft Tissue Tension** – Some shuck is normal when applying a longitudinal distraction force to the hip. Shuck should not be excessive, and the hip should not dislocate (Figure 8A). Rectus femoris tightness (hip in extension, knee flexed) should be no tighter than pre-op.
- 2 Anterior Stability** – Place the leg in full abduction, full extension and hyperextension, while exerting an external rotation force. If the hip cannot be fully extended, it may be too tight. If it dislocates easily, it is too loose and impingement must be addressed or component malposition exists (Figure 8B).
- 3 Posterior Stability** – Place the leg in neutral adduction and 90° flexion. Gradually rotate internally. If it dislocates with minimal internal rotation, it is too loose and impingement must be addressed or component malposition exists (Figure 8C).
- 4 Sleep Position** – Place the leg in the “sleep position” with the operated leg semiflexed, adducted and internally rotated over the other leg. Apply axial force to try to dislocate. This position represents a dangerous unstable position that may be adopted by a patient sleeping on their nonoperated side (Figure 8D).
- 5 Combined Component Positioning** – Place the leg in neutral extension and adduction. Internally rotate the hip 45°. The cup should cover the “northern hemisphere” of the head. This position is an additional test of the positioning of the components in relation to each other.



Figure 8A




Figure 8B



Figure 8C



Figure 8D

Implant insertion

Step 9

Assemble the threaded stem inserter by inserting the stem inserter pommel through the stem inserter frame. Stand the stem inserter upright so that the threaded tip is pointed up (Figure 9A). Screw the implant onto the threaded tip as far as possible.

Flip the assembly over so that the stem tip is now pointing down (Figure 9B). Engage the frame tine into the slots adjacent to the threaded hole on the stem. Screw the pommel until assembly is secure (Figure 9C). Fully tighten the pommel before impaction.

Stem version can be accurately measured by attaching the anteversion handle to the stem inserter frame. This not only allows accurate visualization of anteversion, but will help control rotation of the stem during impaction.



Figure 9A



Figure 9B



Figure 9C

Implant insertion

Step 10

Insert the implant into the canal with hand pressure and verify proper implant version. Use firm mallet blows to seat the implant to the desired level (Figure 10).

Caution Do not use the Stem Inserter Pommel as a stand alone instrument, either for stem insertion or removal

Note Once the implant flutes have engaged the bone, the implant version cannot be changed without removing the implant. The implant can be removed by striking the underside of the threaded stem driver with a mallet.

Surgeon tip When using a porous coated cylindrical stem, some surgeons prefer to know the exact dimensions of the reamers and implants used. In these situations, surgeons may use a ring gauge to measure the reamers and implants to within 0.5mm. Based on those measurements, they may adjust their reaming to tailor implant fit.




Figure 10

Implant trialing and femoral head assembly

Step 11

Once the implant is fully seated, perform a final trial reduction to determine appropriate neck length. Place the desired trial femoral head on the implant and reduce the hip to assess stability and range of motion (Figure 11).



Figure 11

Step 12

Clean and dry the taper with a sterile cloth, place the prosthetic femoral head on the neck taper and firmly impact several times with a femoral head impactor and a mallet (Figure 12).



Figure 12

Conclusion

Irrigate the wound and maintain appropriate hemostasis. If desired, repair the hip capsule and overlying tendons. Close the deep fascial layers, subcutaneous tissue and skin.



Catalog information



Standard offset



High offset

Femoral Implants

Size	Standard offset	High offset	Size	Standard offset	High offset
11	7134-1011	—	16	7134-1016	7134-1026
12	7134-1012	7134-1022	17	7134-1017	7134-1027
13	7134-1013	7134-1023	18*	7134-1018	7134-1028
14	7134-1014	7134-1024	19*	7134-1019	7134-1029
15	7134-1015	7134-1025			

*18 and 19 available upon request



OXINIUM[®] 12/14 Taper Femoral Heads

Neck length	22mm	26mm	28mm	32mm	36mm
-3	—	—	7134-2803	7134-3203	7134-3603
+0	7134-2200	7134-2600	7134-2800	7134-3200	7134-3600
+4	7134-2204	7134-2604	7134-2804	7134-3204	7134-3604
+8	7134-2208	7134-2608	7134-2808	7134-3208	7134-3608
+12	7134-2212	7134-2612	7134-2812	7134-3212	7134-3612
+16	—	—	7134-2816	7134-3216	—

*7134-2340 OXINIUM 40mm Modular Femoral Head

*7134-2344 OXINIUM 44mm Modular Femoral Head



BioloX[™] forte Ceramic Femoral Heads 12/14 Taper

Neck length	28mm	32mm	36mm
+0 (short)	71330280	71330320	71332084
+4 (medium)	71330284	71330324	71332085
+8 (long)	71330288	71330328	71332086



BioloX delta Ceramic Femoral Heads 12/14 Taper

Neck length	28mm	32mm	36mm
+0 (short)	76539160	76539165	71346004
+4 (medium)	76539161	76539166	71346005
+8 (long)	76539162	76539167	71346006

Catalog information



CoCr 12/14 Taper Femoral Heads Cobalt Chromium – ASTM F 799

Neck length	22mm	26mm	28mm	32mm	36mm
-3	—	—	7130-2803	7130-3203	7130-3603
+0	7130-2200	7130-2600	7130-2800	7130-3200	7130-3600
+4	7130-2204	7130-2604	7130-2804	7130-3204	7130-3604
+8	7130-2208	7130-2608	7130-2808	7130-3208	7130-3608
+12	7130-2212	7130-2612	7130-2812	7130-3212	7130-3612
+16	—	—	7130-2816	7130-3216	—

*7134-2640 CoCr 40mm Modular Femoral Head

*7134-2644 CoCr 44mm Modular Femoral Head



Global Femoral Head Trial 12/14 Taper

Length	Neck Color	22mm	28mm	32mm	36mm	40mm	44mm
XS/-3	Green	—	7510-0843	7510-0849	7510-0855	7510-0868	7510-0873
S/+0	Rust	7510-0839	7510-0844	7510-0850	7510-0856	7510-0869	7510-0874
M/+4	Brown	7510-0840	7510-0845	7510-0851	7510-0857	7510-0870	7510-0875
L/+8	Gray	7510-0841	7510-0846	7510-0852	7510-0858	7510-0871	7510-0876
XL/+12	Blue	7510-0842	7510-0847	7510-0853	7510-0859	7510-0872	7510-0877
XXL/+16	Black	—	7510-0848	7510-0854	—	—	—



Titanium Modular Neck Sleeve 12/14 Taper

Neck length		Neck length	
-4	71344245	+4	71344248
+0	71344247	+8	71344249

Use with 40mm and 44mm Oxinium and CoCr Femoral Heads

Box Osteotome
Cat. No. 7136-4002



T-handle
(2 per set)
Cat. No. 7136-4006



Anteversio Handle
(2 per set)
Cat. No. 7136-4012



Osteotomy Guide
Cat. No. 7136-4100



Femoral Canal Finder
Cat. No. 7136-4001



Broach Handle
(2 per set)
Cat. No. 7136-4007



Proximal Reamer
Cat. No. 7136-4015



Catalog information

Trial Neck

Size	Standard offset	High offset
11	7136-7201	—
12	7136-7201	7136-7221
13-14	7136-7202	7136-7222
15-17	7136-7203	7136-7223
18-19	7136-7204	7136-7224



Stem Inserter Pommel

Cat. No. 7136-4011



Stem Inserter Frame

Cat. No. 7136-4008



Broach

Cat. No.	Size	Cat. No.	Size
7136-7010	10	7136-7016	16
7136-7011	11	7136-7017	17
7136-7012	12	7136-7018	18
7136-7013	13	7136-7019	19
7136-7014	14	7136-7020	20
7136-7015	15		



Calcar Reamer

Cat. No. 7136-4004



Femoral Head Impactor


Cat. No. 7136-4009



ECHELON[◇] Primary Instrument Set (Cat. No. ECH002)

ECHELON Starter Tray

Cat. No. 7136-6001





Cat. No.	Description	Ref.
7136-4100	Osteotomy Guide	A
7136-4002	Box Osteotome	B
7136-4001	Femoral Canal Finder	C
7136-4006	T-Handle (2 per set)	D
7136-4015	Proximal Reamer	E
7136-4007	Broach Handle (2 per set)	F
7136-4012	Anteverson Handle (2 per set)	G

ECHELON[◇] Primary Instrument Set (Cat. No. ECH002)

ECHELON Rigid Reamer Tray

Cat. No. 7136-6002




Cat. No.	Description	Ref.
7135-0090	Rigid Reamer, Size 9 mm	A
7135-0095	Rigid Reamer, Size 9.5 mm	B
7135-0100	Rigid Reamer, Size 10 mm	C
7135-0105	Rigid Reamer, Size 10.5 mm	D
7135-0110	Rigid Reamer, Size 11 mm	E
7135-0115	Rigid Reamer, Size 11.5 mm	F
7135-0120	Rigid Reamer, Size 12 mm	G
7135-0125	Rigid Reamer, Size 12.5 mm	H
7135-0130	Rigid Reamer, Size 13 mm	I
7135-0135	Rigid Reamer, Size 13.5 mm	J
7135-0140	Rigid Reamer, Size 14 mm	K
7135-0145	Rigid Reamer, Size 14.5 mm	L

Cat. No.	Description	Ref.
7135-0150	Rigid Reamer, Size 15 mm	M
7135-0155	Rigid Reamer, Size 15.5 mm	N
7135-0160	Rigid Reamer, Size 16 mm	O
7135-0165	Rigid Reamer, Size 16.5 mm	P
7135-0170	Rigid Reamer, Size 17 mm	Q
7135-0175	Rigid Reamer, Size 17.5 mm	R
7135-0180	Rigid Reamer, Size 18 mm	S
7135-0185	Rigid Reamer, Size 18.5 mm	T
7135-0190	Rigid Reamer, Size 19 mm	U
7135-0195	Rigid Reamer, Size 19.5 mm	V
7135-0200	Rigid Reamer, Size 20 mm	W

ECHELON® Broach Tray

Cat. No. 7136-6003



Cat. No.	Description	Ref.
7136-7010	Broach, Size 10	A
7136-7011	Broach, Size 11	B
7136-7012	Broach, Size 12	C
7136-7013	Broach, Size 13	D
7136-7014	Broach, Size 14	E
7136-7015	Broach, Size 15	F
7136-7016	Broach, Size 16	G
7136-7017	Broach, Size 17	H
7136-7018	Broach, Size 18	I
7136-7019	Broach, Size 19	J
7136-7020	Broach, Size 20	K
7136-7103	Revision Trial Neck, Size 18-20	L
7136-7102	Revision Trial Neck, Size 13-17	M
7136-7101	Revision Trial Neck, Size 11-12	N
7136-7201	Trial Neck, Size 10-12, Standard Offset	O
7136-7202	Trial Neck, Size 13-14, Standard Offset	P
7136-7203	Trial Neck, Size 15-17, Standard Offset	Q
7136-7221	Trial Neck, Size 12, High Offset	R
7136-7222	Trial Neck, Size 13-14, High Offset	S
7136-7223	Trial Neck, Size 15-17, High Offset	T
7136-7204	Trial Neck, Size 18-19, Standard Offset	U
7136-7224	Trial Neck, Size 18-19, High Offset	V
7136-4008	Stem Inserter Frame	W
7136-4011	Stem Inserter Pommel	X
7136-4009	Femoral Head Impactor	Y
7136-4004	Calcar Reamer	Z

Total Hip Systems *continued*

CAUTION: Federal law [USA] restricts the subject total hip arthroplasty devices to sale by or on the order of a physician.

Intraoperative *continued*

Common errors include misplacement of the incision, inadequate exposure or mobilization of the femur, inadequate removal of ectopic bone, improper positioning of components, or inadequate proximal support of the femoral component. Studies have indicated a higher risk of implant fatigue fracture in cases with inadequate proximal bone stock or where extended trochanteric osteotomies have been performed. In these cases, it is imperative that adjunctive reinforcement procedures such as bone grafting, cortical strut allografts, cables, and trochanteric plates are utilized to provide adequate proximal support to the femoral component. The use of larger prostheses may also reduce the risk of avoiding prosthetic fatigue fracture. Although these adjunctive reinforcement procedures may minimize the risk of implant failure, they do not ensure a predictable clinical result.

- Prior to closure, the surgical site should be thoroughly cleaned of cement, bone chips, ectopic bone, or other foreign matter. Ectopic bone and/or bone spurs may lead to dislocation and painful or restricted motion.
- Range of motion should be thoroughly assessed for early impingement or joint instability. Postoperative instability (i.e., dislocation) is a leading complication associated with revision surgery and may result in additional surgery.
- Proper positioning of the components is important to minimize impingement which could lead to early failure, premature wear, device related noise, and/or dislocation, all of which may lead to revision surgery.
- To minimize the risk of acetabular shell loosening in uncemented applications, surgeons should consider the use of orthopedic bone fixation devices such as bone screws, spikes, pegs, fins, or other bone fixation devices. To minimize the risk of loose cemented acetabular shells, care should be taken to prevent movement of the implant components while the cement cures.
- Physicians should consider component malposition, component placement, and the effect on range of motion and stability when using modular heads (with sleeves or skirts) and overhang liners.
- For computer assisted surgery systems, it is extremely important to correctly select input parameters (e.g., bony landmarks). Operation of this equipment should be familiar with the anatomy relevant to the procedure. Failure to provide proper input could cause problems such as violation of critical anatomical structures and malpositioned implants, which may lead to revision surgery.
- Trial instrumentation may be provided for the intraoperative assessment of the final implant fit. Do NOT implant trial components.
- Do not implant HA-coated devices in bone cement.
- Inappropriate use of taper sleeves may lead to implant failure which may lead to revision surgery. Select the appropriate sleeve on the Compatible Sleeve Combinations Charts located in the Device Description section of this document.

Postoperative

- Postoperative warnings, precautions, and patient care instructions presented by the physician are extremely important. Gradual weight bearing begins after surgery in ordinary total hip arthroplasty procedures. However, with the trochanter osteotomy or certain complex cases, the weight bearing status should be individualized with the non or partial weight bearing period extended.
- Patients should be warned against unassisted activity, particularly use of toilet facilities and other activities requiring excessive motion of the hip, as they may result in subluxation or dislocation.
- Handle patients with extreme care. Support should be provided to the operative leg when moving the patient. While placing the patient on bedpans, changing dressings and clothing, and similar activities, precautions should be taken to avoid placing an excessive load on the operative leg.
- Postoperative therapy, prescribed by the physician, should be structured to regain muscle strength around the hip and to attain a gradual increase of activities.
- Periodic x-rays, prescribed by the physician, are recommended for comparison to immediate postoperative conditions to detect long-term evidence of changes in position, loosening, bending, and/or cracking of components or loss of bone. If these conditions are evident, patients should be closely observed, the possibilities of further deterioration evaluated, and the benefits of early revision considered.
- If the ceramic head must be revised for any reason and the hip stem is firmly fixed, the revision should be made with a CoCr head and corresponding polyethylene liner and metal shell. If the REFLECTION Ceramic Liner requires revision, both the ceramic liner and the REFLECTION Acetabular Shell cannot be reassembled to any liner. If the R3 poly liner requires revision, and the R3 Acetabular Shell is well-fixed, a new R3 poly liner may be assembled to the existing R3 acetabular shell. If fractured ceramic material is encountered intraoperatively, remove all loose, identifiable fragments and thoroughly irrigate and suction the operative site.
- Prophylactic antibiotics should be recommended to the patient, similar to those suggested by the American Heart Association, for conditions or situations which may result in bacteremia.
- Normal daily activity may be resumed at the physician's direction. Patients should be advised to seek a medical opinion(s) before entering potentially adverse environments that could affect the performance of the implant, such as electromagnetic or magnetic fields, including a magnetic resonance environment.
- The patient should be advised to report any pain, decrease in range of motion, swelling, fever, squeaking, clicking, popping, grating, or grinding noises, and unusual incidences. Patient reports of squeaking, clicking, popping, grating, or grinding should be carefully evaluated as they may indicate position changes in the components which may compromise the durability of the implants.
- Postoperative subluxation may result in higher wear and implant damage.

Cleaning and Sterilization

Cleaning

Refer to the document "Instructions for care, maintenance, cleaning and sterilization of Smith & Nephew orthopaedic devices." This document, reference number 71381339, is available from customer service or via the Smith & Nephew website.

Sterilization

Refer to the product label for the method of sterilization. If not specifically labeled sterile, components are supplied non-sterile and must be cleaned and sterilized prior to surgery. Refer to the document, "Recommendations for decontamination instructions for care, maintenance, cleaning, and sterilization of Smith & Nephew orthopaedic devices," which is available from customer service or via the Smith & Nephew website, for further information regarding the cleaning instructions and the validated sterilization procedures.

Recommended Steam Sterilization Cycle Parameters for Reusable Instruments

- **Dynamic Air Removal (Prevacuum) Steam Cycle:**
 - Exposure temperature: 132°C (270°F); Exposure time: 4 minutes
 - Exposure temperature: 135°C (275°F); Exposure time: 3 minutes
 - Minimum drying time: Wrapped devices - 15 minutes
Containerized devices - 30 minutes
- **Gravity Displacement Steam Cycle:**
 - Exposure temperature: 132°C (270°F)
 - Exposure time: 15 minutes for instruments not in a containment device
30 minutes* for devices in a containment device
 - Minimum drying time: 30 minutes
 - *This sterilization cycle is not considered by the United States Food and Drug Administration (US FDA) to be a standard sterilization cycle. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the selected sterilization cycle specifications (time and temperature).
- **Immediate Use Steam Sterilization or Flash Steam Cycle:**
 - Exposure temperature: 132°C (270°F)
 - Exposure time: Dynamic air removal (pre-vacuum): 4 minutes

For Non-US Customers

- **United Kingdom Steam Cycle:**
 - Pre-vacuum Cycle**
 - Exposure temperature: 134°C (273°F)
 - Exposure time: 3 minutes
 - Vacuum drying time: 30 minutes
- **World Health Organization (WHO) Steam Cycle:**
 - Exposure temperature: 134°C (273°F)
 - Exposure time: 18 minutes
 - Vacuum drying time: 30 minutes

Note: Sterilization evacuation and pulsing should be carried out in accordance with HTM 2010.

Magnetic Resonance Imaging (MRI) Safety

Smith & Nephew hip systems have not been reviewed by the FDA for safety and compatibility in the MR environment. Hip system components have not been tested for heating or migration in the MR environment. Known risks of exposing implant devices to the MR environment include displacement, torque, and radio frequency induced heating. Implant devices may also create image artifacts in MR scans.

Retrieval and Analysis of Harvested Implants

The most important part of surgical implant retrieval is preventing damage that would render scientific examination useless. Special care should be given to protect the implant from damage during handling and shipment. Specifically, for conventional polyethylene or XLPE, use alternative sterilization method other than steam autoclave. Follow internal hospital procedures for the retrieval and analysis of implants harvested during surgery. When handling the harvested implants, use precautions to prevent spread of bloodborne pathogens.

Information

For further information on the medical devices, the information presented herein, or assistance in returning product contact Customer Service at (800) 238-7538 for calls within the continental USA and (901) 396-2121 for all international calls.

Smith & Nephew Inc.
1450 Brooks Road
Memphis, TN 38116 U.S.A.
Tel.: 901-396-2121
www.smith-nephew.com

Smith & Nephew Orthopaedics GmbH
Allemannenstrasse 14
78532 Tuttlingen, Germany
Tel.: 07462/208-0
Fax: 07462/208-155

Explanation of symbols used in labeling:

- H₂O₂ – Hydrogen peroxide sterilization
- ID – Inner diameter
- OD – Outer diameter
- S/+0 – Short
- M/+4 – Medium
- L/+8 – Long
- SO – Standard offset
- H or HO – High offset



– For use with bone cement



– For use without bone cement

*Trademark of Smith & Nephew. Certain marks registered U.S. Patent and Trademark Office. All trademarks acknowledged.



Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA



Smith & Nephew Orthopaedics GmbH
Allemannenstrasse 14
78532 Tuttlingen, Germany



www.smith-nephew.com
+1 800 238 7538 U.S. Customer Service
+1 901 396 2121 International Customer Service

Smith & Nephew, Inc.
1450 Brooks Road
Memphis, TN 38116
USA

www.smith-nephew.com

Telephone: 1-901-396-2121
Information: 1-800-821-5700
Orders and Inquiries: 1-800-238-7538